

Table Drape Scatter Shields

Perfect for Cath Labs, Interventional Suites, and any imaging office that has a C-Arm. These single sided drapes do not interfere with c-arm placement and are perfect for lowering those high dosimeter badge readings by lowering scatter radiation to healthcare staff. The pleated panels overlap at least 1" to provide that maximum coverage you can depend on. Offering several common sizes to choose from.

Manufactured using 0.5mm Pb equivalent lightweight lead instead of the older regular lead making it much easier to use and attach to the imaging table. Our 2" non-media table flap has a sewn in reinforced "hook" that will attach on top of a 2" wide self-adhesive black "loop" that can be secured to a table or gurney. These can easily be cleaned with 10% bleach solution or any medical grade disinfectant wipe or spray.

For facilities needing double sided protection, simply order two drapes (one for each side of the table/gurney). Anesthesia drapes are also available.

- **Core Material: Lightweight Lead**
- **Core Material Thickness: 0.5mm Pb equivalent protection**
- **Covering: Classic Easy to Clean Nylon**
- **Colors: Black**
- **Three standard widths: 24", 30", 36" and 48"**
- **Three standard heights: 24" , 30" , 36" , and 42"**
- **Manufactured in the USA**
- **10-15 business days for manufacturing**

Custom table drapes are also available for tables needing "tab loops" instead of self adhesive loop to connect drapes to tables with bars on the side. Custom table drape sizes (length and width) are also available.



Table Drape Scatter Shields

Directions for Attachment to Your Imaging Tables

4 Easy Steps...



Step 1: Remove the clear plastic from the self-adhesive “loop” strip that is supplied with your table drape.



Step 2: Firmly press the self-adhesive “loop” strip onto the edge of your table top.



Step 3: Carefully keep removing the clear plastic backing from the self adhesive strip.



Step 4: Attach your table drape to the self-adhesive “loop” and firmly press down to secure the drape in place.



Table Drape Scatter Shields Cleaning & Care Instructions



Which Products Should You Use to Clean Your Drapes?

Many cleaners deteriorate the color and the core material. Over time, this may also cause tearing or holes, which compromises the safety of the table scatter drape. These cleaning products focus on disinfection while maintaining the integrity of the garment. These cleaners will not harm or deteriorate your drapes and have been tested thoroughly with many applications.

Approved Cleaners:

- [Revolution Scrubbles Lead Apron & Table Drape Cleaner](#)
- [Clorox Hydrogen Peroxide Cleaner Disinfectant Spray](#)
- [Clorox Hydrogen Peroxide Cleaner Disinfectant Wipes](#)



How to Clean Your Table Drape

With [Clorox Hydrogen Peroxide Cleaner Disinfectant Wipes](#), start by scrubbing the table drape panels in a circular motion. This helps ensure you are covering the full surface area where bacteria may live. Continue this process until you have wiped down the entire table drape and each panel.

When using the [Revolution Scrubbles Lead Apron and Table Drape Cleaner](#), or the [Clorox Hydrogen Peroxide Cleaner Disinfectant Spray](#), spray the disinfectant and scrub thoroughly in circular motions. Use the brush to give your table drape a much deeper clean. Continue this process until you have scrubbed the entire drape.



How to Store Your Table Drape



Your table drape can be left on the imaging table or you can purchase a table drape storage bag to store your table drape when it is not in use.
Sku# **GDN-DUFFLE-ZIP**

